



# Proteogenix

**Proteogenix SAS**

www.proteogenix.science  
19 rue de la Haye  
67300 Schiltigheim, France  
Tel +33 (0) 390 20 54 70  
Fax +33 (0) 978 53 36 90  
products@proteogenix.fr

## Batiraxcept Biosimilar - Anti-AXLLG mAb - Research Grade

### OVERVIEW

---

**PRODUCT REFERENCE** PX-TA2004

**CLONALITY**

**ISOTYPE** Fusion - [AXL (AXL receptor tyrosine kinase, ARK, JTK11, Tyro7, UFO)]<sub>2</sub> - IGHG1 Fc  
(Fragment constant)

**SOURCE**

### PRODUCT INFORMATION

---

**ORIGIN SPECIES**

**EXPRESSION SYSTEM** XtenCHO

**BUFFER** 0.01M PBS, pH 7.4.

**STORAGE CONDITIONS** Store at 4 °C for short term. Freeze and store at -80 °C for long term.  
Addition of 20 to 40 % glycerol may improve cryoprotection.  
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.  
Aliquot the product according to your needs to avoid freeze/thaw cycles.

**ALIAS/SYNONYMS** anti-AXLLG, Growth arrest-specific protein 6, AXL receptor tyrosine kinase ligand,  
GAS6, GAS-6



**WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.**