



Nexviazyme/Avalglucosidase alfa Biosimilar - Anti-GAA mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA2028
CLONALITY	Monoclonal Antibody
ISOTYPE	0
SOURCE	

PRODUCT INFORMATION

ORIGIN SPECIES

EXPRESSION SYSTEM XtenCHO

BUFFER 0.01M PBS, pH 7.4.

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.

ALIAS/SYNONYMS anti-GAA, Lysosomal alpha-glucosidase, Acid maltase, Aglucosidase alfa

