



# Proteogenix

**Proteogenix SAS**

www.proteogenix.science  
19 rue de la Haye  
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70  
Fax +33 (0) 978 53 36 90  
products@proteogenix.fr

## Albiglutide Biosimilar - Anti-GLP-1 receptor mAb - Research Grade

### OVERVIEW

---

**PRODUCT REFERENCE** PX-TA2029

**CLONALITY**

**ISOTYPE** Fusion - (GCG (glucagon) 98-127 (glucagon-like peptide 1, GLP1) 7-36)<sub>2</sub> - ALB  
(albumin, human serum albumin, HSA) 25-609

**SOURCE**

### PRODUCT INFORMATION

---

**ORIGIN SPECIES** Homo sapiens

**EXPRESSION SYSTEM** XtenCHO

**BUFFER** 0.01M PBS, pH 7.4.

**STORAGE CONDITIONS** Store at 4 °C for short term. Freeze and store at -80 °C for long term.  
Addition of 20 to 40 % glycerol may improve cryoprotection.  
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.  
Aliquot the product according to your needs to avoid freeze/thaw cycles.

**ALIAS/SYNONYMS** anti-GLP-1 receptor, Glucagon-like peptide 1 receptor, GLP1R, GLP-1-R, GLP-1R



**WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.**