



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Albinterferon alfa-2b Biosimilar - Anti-hepatitis C virus fusion protein - Research Grade

OVERVIEW

PRODUCT REFERENCE PX-TA2030

CLONALITY

ISOTYPE Fusion - [ALB (albumin, human serum albumin, HSA) 25-609 - IFNA2 (interferon alpha 2) *b (R46,H57) 24-188]

SOURCE

PRODUCT INFORMATION

ORIGIN SPECIES Homo sapiens

EXPRESSION SYSTEM XtenCHO

BUFFER 0.01M PBS, pH 7.4.

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.

ALIAS/SYNONYMS anti-hepatitis C virus



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.