



Albutrepenonacog alfa Biosimilar - Anti-hemophilia fusion protein - Research Grade

OVERVIEW

PRODUCT REFERENCE PX-TA2031

CLONALITY

ISOTYPE Fusion - [F9 (coagulation factor 9, coagulation factor IX) recombinant - ALB (albumin, human serum albumin, HSA) 25-609]

SOURCE

PRODUCT INFORMATION

ORIGIN SPECIES Homo sapiens

EXPRESSION SYSTEM XtenCHO

BUFFER 0.01M PBS, pH 7.4.

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.

ALIAS/SYNONYMS anti-Plasma thromboplastin component, F9, Christmas factor, Coagulation factor IX, PTC

