



Asfotase alfa Biosimilar - Anti-TNSALP mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE PX-TA2032

CLONALITY

ISOTYPE Fusion - [ALPL (alkaline phosphatase, TNSALP, tissue non-specific alkaline phosphatase) 18-502]2 - IGHG1 Fc (Fragment constant) - Asp-Ile-deca-aspartate

SOURCE

PRODUCT INFORMATION

ORIGIN SPECIES Homo sapiens

EXPRESSION SYSTEM XtenCHO

BUFFER 0.01M PBS, pH 7.4.

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.

ALIAS/SYNONYMS anti-TNSALP, Alkaline phosphatase liver/bone/kidney isozyme, AP-TNAP, Alkaline phosphatase, tissue-nonspecific isozyme, ALPL

