

Anzurstobart Biosimilar - Anti-BIT mAb - Research Grade

OVERVIEW

| | |
|--------------------------|---------------------|
| PRODUCT REFERENCE | PX-TA2044 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | IgG1-kappa |
| SOURCE | |

PRODUCT INFORMATION

| | |
|---------------------------|---|
| ORIGIN SPECIES | Homo sapiens |
| EXPRESSION SYSTEM | XtenCHO |
| BUFFER | 0.01M PBS, pH 7.4. |
| STORAGE CONDITIONS | Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles. |
| ALIAS/SYNONYMS | anti-BIT, Signal-regulatory protein alpha-1, Brain Ig-like molecule with tyrosine-based activation motifs, MFR, SIRP, Sirp-alpha-1, Sirp-alpha-2, PTPNS1, SIRPA, SHP substrate 1, CD172 antigen-like family member A, Sirp-alpha-3, Inhibitory receptor SHPS-1, |

