



Ralzapastotug Biosimilar - Anti-VSTM3 mAb - Research Grade

OVERVIEW

| | |
|--------------------------|---------------------|
| PRODUCT REFERENCE | PX-TA2051 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | IgG1-kappa |
| SOURCE | |

PRODUCT INFORMATION

| | |
|---------------------------|---|
| ORIGIN SPECIES | Homo sapiens |
| EXPRESSION SYSTEM | XtenCHO |
| BUFFER | 0.01M PBS, pH 7.4. |
| STORAGE CONDITIONS | Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles. |
| ALIAS/SYNONYMS | anti-VSTM3, TIGIT, V-set and transmembrane domain-containing protein 3, VSIG9, V-set and immunoglobulin domain-containing protein 9, T-cell immunoreceptor with Ig and ITIM domains |

