



Prafnosbart Biosimilar - Anti-SKR1 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA2084
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa
SOURCE	

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	XtenCHO
BUFFER	0.01M PBS, pH 7.4.
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-SKR1, Activin receptor type-1, Activin receptor-like kinase 2, TSR-I, ACVRLK2, ACTR-I, ALK-2, Serine/threonine-protein kinase receptor R1, TGF-B superfamily receptor type I, ACVR1, Activin receptor type I

