



Tulisokibart Biosimilar - Anti-TNFSF15 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA2090
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa
SOURCE	

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	XtenCHO
BUFFER	0.01M PBS, pH 7.4.
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	anti-TNFSF15, Vascular endothelial cell growth inhibitor, VEGI, TL1, TNF ligand-related molecule 1, Tumor necrosis factor ligand superfamily member 15

