



ProteoGenix

ProteoGenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

ProteoGenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME Onvapegleukin Alfa Biosimilar - Anti-IL2 mAb - Research Grade

PRODUCT REFERENCE PX-TA2158

CLONALITY

ISOTYPE Human interleukin-2, variant (R38>C, C125>S)

SOURCE CAS: 2764760-02-1

PRODUCT INFORMATION

ORIGIN SPECIES

PRODUCTION HOST XtenCHO

FORM Liquid

BUFFER 0.01M PBS, pH 7.4.

STABILITY & STORAGE Store at 4°C for short term. Freeze and store at -20°C or -80°C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Avoid freeze/thaw cycles, aliquot the product according to your needs.

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY