

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Brivekimig Biosimilar - Anti-CD252,TNFA,ALB mAb - Research Grade
PRODUCT REFERENCE	PX-TA2180
CLONALITY	Monoclonal Antibody
ISOTYPE	Ig IG single chain of 5 VH (VH-VH-VH'-VH''-VH') monomer
SOURCE	Trispecific, Pentavalent, CAS: 2834733-45-6

PRODUCT INFORMATION

ORIGIN SPECIES	Human
PRODUCTION HOST	XtenCHO
FORM	Liquid
BUFFER	0.01M PBS, pH 7.4
STABILITY & STORAGE	Store at 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.