

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Imelciment Biosimilar - Anti-Vixticibart mAb - Research Grade
PRODUCT REFERENCE	PX-TA2195
CLONALITY	Monoclonal Antibody
ISOTYPE	Ig (H-gamma4_L-kappa)_G4hCH2CH3
SOURCE	Idiotypic, CAS: 2845128-03-0

PRODUCT INFORMATION

ORIGIN SPECIES	Human
PRODUCTION HOST	XtenCHO
FORM	Liquid
BUFFER	0.01M PBS, pH 7.4
STABILITY & STORAGE	Store at 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.