

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Pegrizeprium Biosimilar - Anti-CD28 mAb - Research Grade
<b>PRODUCT REFERENCE</b>	PX-TA2210
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	Ig Fab VH-G1CH1h_L-kappa
<b>SOURCE</b>	CAS: 2829316-24-5

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	Human
<b>PRODUCTION HOST</b>	XtenCHO
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS, pH 7.4
<b>STABILITY &amp; STORAGE</b>	Store at 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.