



# ProteoGenix

**ProteoGenix S.A.S**

19 rue de la Haye  
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

[products@proteogenix.fr](mailto:products@proteogenix.fr)

**ProteoGenix INC.**

1000 Bearcat Way, Suite  
100, Morrisville, NC 27560,  
USA

Tel +1(919)234-1277

[products@proteogenix.fr](mailto:products@proteogenix.fr)

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Podentamig Biosimilar - Anti-TNFRSF17,ALB,CD3e mAb - Research Grade
<b>PRODUCT REFERENCE</b>	PX-TA2211
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	Ig single chain VH-VH'-scFvhl
<b>SOURCE</b>	Trispecific, CAS: 2851863-52-8

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	Human
<b>PRODUCTION HOST</b>	XtenCHO
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS, pH 7.4
<b>STABILITY &amp; STORAGE</b>	Store at 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

You can contact us at [products@proteogenix.fr](mailto:products@proteogenix.fr) for any question

**OUR PRODUCTS ARE FOR RESEARCH USE ONLY**