



ProteoGenix

ProteoGenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

ProteoGenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Invikafusp Biosimilar - Anti-IL2 fusion protein - Research Grade
PRODUCT REFERENCE	PX-TA2230
CLONALITY	
ISOTYPE	Humanized IgG1-kappa anti-human TCR V α 6/V β 10 monomer disulfide bridged to human IL-2 fragment fused via (G4S) ₂ peptide linker to Ig G1 Fc fragment.
SOURCE	CAS: 2894138-60-2

PRODUCT INFORMATION

ORIGIN SPECIES	Human
PRODUCTION HOST	XtenCHO
FORM	Liquid
BUFFER	0.01M PBS, pH 7.4
STABILITY & STORAGE	Store at 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY