



Proteogenix

Proteogenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

Proteogenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|--|
| PRODUCT NAME | Vatreptacog alfa Biosimilar - Anti-FVIIa Mutant Version (FVIIaDVQ) fusion protein - Research Grade |
| PRODUCT REFERENCE | PX-TA2274 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | A recombinant factor VIIa (rFVIIa) analog with three amino acid substitutions (Asp158Val, Val296Glu, and Glu298Met). |
| SOURCE | P08709 |

PRODUCT INFORMATION

ORIGIN SPECIES

PRODUCTION HOST

FORM Liquid

BUFFER

STABILITY & STORAGE Store at 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY