



ProteoGenix

ProteoGenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90

products@proteogenix.fr

ProteoGenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|---|
| PRODUCT NAME | Mirostipen Biosimilar - Anti-CD191 - Research Grade |
| PRODUCT REFERENCE | PX-TA2279 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | CCL23 protein |
| SOURCE | P32246 |

PRODUCT INFORMATION

ORIGIN SPECIES

PRODUCTION HOST

FORM Liquid

BUFFER

STABILITY & STORAGE Store at 4°C for short term. Freeze and store at -20°C or -80°C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Avoid freeze/thaw cycles, aliquot the product according to your needs.

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY