

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Sozinibercept Biosimilar - Anti-VEGF-C fusion protein - Research Grade
PRODUCT REFERENCE	PX-TA2286
CLONALITY	Monoclonal Antibody
ISOTYPE	Vascular Endothelial Growth Factor Receptor 3 fusion protein with IgG1 Fc Fragment.
SOURCE	P49767 & O43915

PRODUCT INFORMATION

ORIGIN SPECIES**PRODUCTION HOST**

FORM Liquid

BUFFER

STABILITY & STORAGE Store at 4°C for short term. Freeze and store at -20°C or -80°C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Avoid freeze/thaw cycles, aliquot the product according to your needs.