



# ProteoGenix

**ProteoGenix S.A.S**

19 rue de la Haye  
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

[products@proteogenix.fr](mailto:products@proteogenix.fr)

**ProteoGenix INC.**

1000 Bearcat Way, Suite  
100, Morrisville, NC 27560,  
USA

Tel +1(919)234-1277

[products@proteogenix.fr](mailto:products@proteogenix.fr)

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	GAA/Agglucosidase Alfa Biosimilar - Anti-GAA - Research Grade
<b>PRODUCT REFERENCE</b>	PX-TA2290
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	A recombinant GAA/Agglucosidase alfa
<b>SOURCE</b>	P10253

### PRODUCT INFORMATION

---

#### ORIGIN SPECIES

#### PRODUCTION HOST

**FORM** Liquid

#### BUFFER

**STABILITY & STORAGE** Store at 4°C for short term. Freeze and store at -20°C or -80°C for long term.  
Addition of 20 to 40 % glycerol may improve cryoprotection.  
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.  
Avoid freeze/thaw cycles, aliquot the product according to your needs.

You can contact us at [products@proteogenix.fr](mailto:products@proteogenix.fr) for any question

**OUR PRODUCTS ARE FOR RESEARCH USE ONLY**