



Proteogenix

Proteogenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

Proteogenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Teriparatide Biosimilar - Anti-PTH - Research Grade
PRODUCT REFERENCE	PX-TA2297
CLONALITY	Polyclonal Antibody
ISOTYPE	Biosynthetic Parathyroid hormone (PTH) analogs.
SOURCE	P01270

PRODUCT INFORMATION

ORIGIN SPECIES

PRODUCTION HOST

FORM Liquid

BUFFER

STABILITY & STORAGE Store at 4°C for short term. Freeze and store at -20°C or -80°C for long term.
Addition of 20 to 40 % glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Avoid freeze/thaw cycles, aliquot the product according to your needs.

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY