

Human CD102/ICAM2 recombinant protein

OVERVIEW

CATALOG NO	PX-P4034
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Partial
MOLECULAR WEIGHT	25.62kDa
DELIVERED WITH TAG	Yes
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
PURITY	80%
UNIPROT ID	P13598
FORM	Liquid
BUFFER	50 mM Tris-HCl pH 8, 150 mM NaCl 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
STABILITY & STORAGE	

More info

SYNONYMS

ICAM2, CD102, CD102 antigen, ICAM 2, Intercellular adhesion molecule 2, Ly60

UNIPROT LINK

<http://www.uniprot.org/uni/prot/P13598>

SEQUENCE

MSSFGRYRTLTVLFTLICCPGSDEKVFVHVRPKKLAVEPKGSLEVNCSTTCNQPEVGGLETSLDKILLDEQAQWKHYLVSNISHDTVLCQ
HFTCSGKQESMNSNVSVYQPPRQVILTLQPTLVAVGKSFTIECRVPTVEPLDSLTLFLFRGNETLHYETFGKAAPAPQEATATFNSTADRED
GHRNFSCLAVLDLMSRGGNIFHKHSAPKMLEIYEPVSDSQGSHHHHHH