



HBV S/L-HBsAg/L glycoprotein Protein, C-S Tag, recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P5456
ORIGIN SPECIES	Hepatitis B virus
FRAGMENT	Met175-Ile400(Pro223Leu, Ala351Thr)
TAG	C-terminal S Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	29.41 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9E6S4
ALIAS/SYNONYMS	

