



CD365 / KIM-1 / TIM-1 / HAVCR1, N-His, recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P5622
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Ser21-Gly295
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	31.80 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q96D42
ALIAS/SYNONYMS	

