



IL6, C-His & amp & Avi tag, recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P5794
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Val30-Met212
TAG	C-terminal His & amp & Avi Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	16.75 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P05231
ALIAS/SYNONYMS	

