



TNF α / TNF-alpha, N-His, recombinant protein (E.coli)

OVERVIEW

PRODUCT REFERENCE	PX-P5960
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Val77-Leu233
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	19.66 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P01375
ALIAS/SYNONYMS	

