



Human ROR1 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6002
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Glu403
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	43.4 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q01973
ALIAS/SYNONYMS	Inactive tyrosine-protein kinase transmembrane receptor ROR1, Neurotrophic tyrosine kinase, receptor-related 1, NTRKR1

