



Human GITR recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6015
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Gln26-Glu161
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	16.84 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9Y5U5
ALIAS/SYNONYMS	Activation-inducible TNFR family receptor, GITR, TNFRSF18, CD357, Tumor necrosis factor receptor superfamily member 18, AITR, Glucocorticoid-induced TNFR-related protein

