

Human GITRL recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6016
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Glu74-Ser199
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Baculovirus-Insect cells
MOLECULAR WEIGHT	13.97 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9UNG2
ALIAS/SYNONYMS	AITRL, GITRL, TL6

