



Human CD328/SIGLEC7 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6023
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Gly18-Leu353
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	39.31 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9Y286
ALIAS/SYNONYMS	QA79 membrane protein, Adhesion inhibitory receptor molecule 1, D-siglec, CDw328, Siglec-7, p75, SIGLEC7, AIRM-1, Sialic acid-binding Ig-like lectin 7, CD328, AIRM1

