



## Human GREM1 recombinant protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P6037
<b>ORIGIN SPECIES</b>	Human
<b>FRAGMENT</b>	Lys26-Asp184
<b>TAG</b>	N-terminal His Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Prokaryotic expression
<b>HOST</b>	Escherichia coli (E.coli)
<b>MOLECULAR WEIGHT</b>	20.46 kDa
<b>FORM</b>	
<b>BUFFER</b>	PBS, pH7.5
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	O60565
<b>ALIAS/SYNONYMS</b>	Cysteine knot superfamily 1, BMP antagonist 1, CKTSF1B1, Down-regulated in Mos-transformed cells protein, Increased in high glucose protein 2, IHG-2, DRM, GREM1, Gremlin-1, Cell proliferation-inducing gene 2 protein, DAN domain family member 2, DAND2

