



Human CSF3 / G-CSF recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6039
ORIGIN SPECIES	Human
FRAGMENT	Thr31-Pro207
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	21.28 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P09919
ALIAS/SYNONYMS	C17orf33, Lenograstim, CSF3, Pluripoietin, Filgrastim, Granulocyte colony-stimulating factor, GCSF, G-CSF

