



Human Flt3 ligand / FLT3LG recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6044
ORIGIN SPECIES	Human
FRAGMENT	Ser25-Pro184
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	21.05 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P49771
ALIAS/SYNONYMS	SL cytokine, FLT3LG, Fms-related tyrosine kinase 3 ligand, Flt3 ligand, Flt3L

