



Human IL17C recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6051
ORIGIN SPECIES	Human
FRAGMENT	His19-Val197
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	22.04 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9P0M4
ALIAS/SYNONYMS	Cytokine CX2, IL17C, IL-17C, Interleukin-17C

